reuse and environment damages.

DO NOT INGEST.

CONTRAINDICATION:

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40x3x3g syringes **ATTAQUE GEL** + accessories.

Reg. ANVISA: 10298550087.

COMPOSITION:

37% Orthophosphoric Acid; Methylparaben; Blue Dyeing (C152015); Thickening and Deionized Water.

INDICATION:

For etching enamel and dentin, allowing better retention of composite material used in restoration, attachment, cementing and sealing.

TECHNICAL INFORMATION:

The 37% Orthophosphoric Acid used for 15 to 30 seconds, is the most efficient and used way to attack dentin or enamel. The use in dentin promotes a fast and total smear layer removal, opening the extremities of the tubules, allowing the resin to penetrate easily; besides, the fact of calcium and phosphate ions, present in dentin and in dentine fluids, are able to react with ortho phosphoric acid, neutralizing its toxic effect on dentin. Therefore, it promotes an increase of micro-mechanical retention and decreases the marginal micro-leakage. The coloring of **ATTAQUE GEL** allows a perfect visualization where to apply it, avoiding any risk of conditioning in other areas. Its viscosity allows a precise placement of the gel.

USE INSTRUCTIONS:

1. After prophylaxis with pumice powder, wash and dry the preparation; Isolation with a rubber dam is recommended;
 2. Apply **ATTAQUE GEL** in the areas to be treated; Time of acting : 15 seconds for Dentin and 30 seconds for enamel;
- Attention: The deep cavities must be protected with calcium hydroxide or glass ionomer cement before the use of **ATTAQUE GEL**;
3. Fluorated tooth can need a second application;
 3. Wash the conditioned area with plenty of water for 30 seconds;
 4. Dry, blowing air (oil free).

SPECIAL CARES:

- Avoid contact with the skin or mucosa. In case of contact, rinse thoroughly with water.
- After the use, withdraw the syringe pressure by pulling back the plunger.
- During product manipulation, the professional should use gloves and safety glasses;
- Protect the patient with apron and safety glasses.
- In case of contact with the eyes, wash with water, and if necessary, search for medical attention
- Protect the neighbor teeth, with polyester strip, in places where you do not want etch.
- Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its