3 application tips.

INDICATION:

MASTER CEMENT LC is a light cured adhesive cementation system, indicated for adhesive cementation of indirect restorations of ceramic, composites or ceromer, inlays, onlays and veneer, up to 2,0mm to the dental structure.

Reg. ANVISA 10298550109

COMPOSITION:

MASTER CEMENT LC (Light Cured): Bisphenol A glicidilmethecrylate; Dimethacrylate groups; Inorganic fillers; Catalyst and Pigments.

TECHNICAL INFORMATIONS:

MASTER CEMENT LC is a light cured adhesive cementation system, indicated to join adhesively indirect restorations up to 2,00mm, to the dental structure. It is compatible with adhesive agents of multiple uses (a.e. Master Bond Multi-Use) as well as with single components adhesives (a.e. Master Bond). Your ceramic fillers attributes to the product great physical qualities. It's presented in different colors, which confers excellent esthetic propriety at veneers and prosthesis without metal cementation.

MASTER CEMENT TRY IN, which may come with the product, is a proof cementation material, based on glycerin and soluble in water. It acts promoting a color simulation which is going to be obtained through the use of **MASTER CEMENT LC**, when completely polymerized. It is presented on the same colors of the resinous cement **MASTER CEMENT LC**.

USE INSTRUCTIONS:

1. Remove of the temporary element and the clean the dental structure with pumice stone and water;

2. Prove the restoration using **MASTER CEMENT TRY-IN**, to a better esthetic result;

Obs.: To avoid the risk of fracture of the restoration, do not check the occlusion in this stage;

3. Remove completely the **MASTER CEMENT TRY-IN** from the restoration with water spray and dry with air jet immune of water and oil;

4. Do the restoration treatment in accordance to the recommended technique of the material in use. Porcelain and composites restorations should be treated in your internal face with **PORCELAIN CONDITIONER** for 1 minute. Then, wash and dry with air jet without oil, and apply **SILANE**. Wait for 30 to 60 seconds and apply a new layer;

5. Promote the etching of the dental structure with **ATTAQUE GEL** (37% Orthophosphoric Acid) and leave acting for 15 seconds on the dentin and for 30 seconds on the enamel. Then, wash with plenty of water and dry without dehydrate the dental surface; OBS.: A light glow of humidity should be observed on the surface.

6. Apply the elected adhesive agent, according to manufacturer's instructions;

OBS.: **MASTER CEMENT LC** is compatible with multiple uses adhesive agents (a.e. **MASTER BOND MULTI-USE**) as well as with single components union agent (a.e. **MASTER BOND**).

7. Apply **MASTER CEMENT LC** inside the restoration with the help of a brush or a spatula, place the restoration in position, keeping it slightly pressed and remove the excess with the help of a proper instrument (a.e. spatula or brush);

8. Light cure for 20 seconds in each segment, starting by the proximal margins;

9. Remove the excess of the light cured material using specific points and do the finishing with silicone points and finishing disks.

SPECIAL CARES:

Eliminate completely all cementation proof material (**MASTER CEMENT TRY IN**) before do the definitive cementation. Residues of this material can damage the adhesion between the restoration and the resinous cement.

Protect the pulp with a lining material or cavitary base if necessary.

Do not use materials containing Eugenol, because Eugenol influences negatively at the polymerization time.

Observe that exists a slightly difference between the shades of **MASTER CEMENT TRY IN** and the correspondent shades of **MASTER CEMENT LC**, before the photo activation, due to the fact that the coloration is changed after the light curing. The shade of **MASTER CEMENT TRY IN** is correspondent to the shade of **MASTER CEMENT LC**, when completely light cured.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRA-INDICATIONS:

MASTER CEMENT LC is contra-indicated for people with related sensitivity to any components of the product.


WARNINGS:

Do not use material after expiry date.

The manufacture is not responsible for any damage resulting from a misuse or not provided for the instructions of use.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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