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IONOFAST

Light cured glass ionomer cement ready to use

biodinâmica®

CONTENTS:

1x5g Syringe of **IONOFAST**;

1x2,5g Syringe of **IONOFAST**;

2x2,5g Syringe of **IONOFAST**;

application tips.

ANVISA Reg. number: 10298559031.

COMPOSITION:

Glass of calcium fluorosilicate, barium and aluminum; Methacrylic monomers; Silicon dioxide; catalyst; stabilizer; pigments.

INDICATION:

IONOFAST is indicated for cavity liner and base under all kind of restorative material (amalgam, ceramic or composites). Pits and fissures sealing and restoration of small cavities. Filing small retention or defects in preparation.

TECHNICAL INFORMATION:

IONOFAST is a monocomponent material, ready to use, using the lightcuring mechanisms for setting; therefore is considered a lightcuring material. The setting reaction of the material occurs by the polymerization of methacrylates and dimethacrylates groups. The presence of methacrylates in the ionomer does not interfere in the interaction between the tooth structure and cement. The adhesion through the ion exchange layer and the fluoride release stays unchanged. The resin component, upon activation by light, protects the cement from early water loss. Thus, the photoactive cements are immediately resistant to water incorporation or water loss and need not sealing. Glass ionomer releases fluoride and maintains the release for a long time, due to the fact that fluoride is not part of the cement matrix constitution. Cement can incorporate fluoride ions available in the oral cavity, provided by topical applications for professional or fluoridated toothpastes for example.

USE INSTRUCTIONS:

1. Prepare the cavity, following the restorative technique. Wash and dry the surface;
2. Apply a cavity base in deep areas (Calcium hydroxide cement);
3. Apply **IONOFAST** in the cavity, directly from the syringe (avoid cross-contamination, decontaminating the syringe tip after use). When a thicker layer is necessary (> 1mm), apply and polymerize each layer.
4. Polymerize at a distance of 2mm for 20 seconds, when using a light curing with a wavelength of 400 to 500nm with a maximum emission peak at 468nm (blue light) and minimum power of 400 mW/cm², total

supply of at least 16J of energy, when the distance of the light curing is less than 5mm;
5. Finish the surface of the ionomer with rotating instruments and restore the cavity according to the technique selected and the instructions of the manufacturer.

SPECIAL CARES:

- After use, keep flasks tightly closed to avoid exposure to light and contaminants;
- Do not storage the product near other products containing Eugenol. It might interfere on lightcuring;
- Eyes contact: rinse immediately with water; seek medical advice if the irritation persists.
- Skin contact: rinse the affected area with water and soap. Consult a doctor if the irritation persists.

CONTRAINDICATION:


IONOFAST should not be used in patients is known to have sensibility to any of the ingredients listed.

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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