USE INSTRUCTIONS:

1. Select the gutta percha point;
2. Put on a flat glass plate 2 or 3 drops of **Liquid PULP FILL** and go adding and spatulating **Powder PULP FILL** during 3 minutes until a thick consistency mass is obtained; The ideal consistency is when moving away the spatula from the mass occurs the formation of a 2,5cm length filament.
3. After gutta percha point is totally involved by the cement, obturate the root canal according with the indicated endodontic technique.

Note: **PULP FILL** can be inserted into the conduct using lentulo drills.

SPECIAL CARES:

It is recommended to avoid **PULP FILL** overflow in the periapical tissue although this product is well tolerated when overflowed.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:


PULP FILL should not be used in patients is known to have sensibility to any of the ingredients listed.

WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.
"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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PULP FILL Liquid: Eugenol and Oil of Sweet Almond.

INDICATION:

PULP FILL is an Endodontic cement to obturate the root canal.

TECHNICAL INFORMATION:

PULP FILL is an endodontic cement, used as intra canal filling material.

The liquid and the powder provide one viscose and consistent mixture, which is used as gutta-percha point. The setting of this cement is caused by the reaction between the Zinc Oxide and Eugenol, which in the presence of water produces the zinc eugenolate. The resin improves the characteristics of the mixture and delays the setting time. The radiopacity is optimized by the bismuth and barium salts.

The formulation of **PULP FILL** is similar to the one developed by GROSSMAN. **PULP FILL** was developed seeking to achieve the requests fixed by GROSSMAN, a perfect endodontic cement: (1) adhesiveness: should be sticky when mixed, to offer a good adhesion between the point and the canal; (2) perfect sealing: should provide a perfect close; (3) radiopaque: should be radiopaque in order to be viewed in the x-ray; (4) small particles: the powder particles should be very small to be mixes easily with the liquid; (5) not suffer contraction: should not contract after insertion into the duct; (6) should not stain the dental structure; (7) antibacterial action: should be bacteriostatic or at least, not allow the bacterial growth; (8) should hard slowly; (9) not soluble: should be insoluble in the liquids; (10) biocompatible: should be tolerated by the tissues; (11) possible to be removed: Should be insoluble in common solvents, in case that is necessary to remove the filling material.