Master Flow

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Master Flow

Light Cured Low Viscosity Composite

CONTENTS

Introductory Kit:
3x2g syringes MASTER FLOW - shades: A2 - A3,5 - OA2.
6 application tips.
Reposition:
1x2g syringe MASTER FLOW - shades: A1 - A2 - A3 - A3,5 - B1 - B2 - B3 - C3 - OA2.
1 application tip.
ANVISA nº 10298550092.

COMPOSITION:

Matrix of Methacrylic Monomers, Inorganic Load, Initiators, Stabilizers and Pigments.

INDICATION:

MASTER FLOW is a composite indicate for anterior or posterior definitive esthetic dental restoration as a first application (initial layer) under light cured composite (a.e. Master Fill); for micro-cavities restoration; for cervical restoration; for preventive restoration in premolar and molar teeth; for ceramic and composites restoration repairs; for dental structures fillings; for the adhesive cementing of ceramic or composite restoration; for splinting of teeth with mobility; for dental fragments bonding.

TECHNICAL INFORMATION:

MASTER FLOW is a light cured, hybrid composite resin with organically modified glass-ceramic microparticles, giving the product greater resistance, greater shine, and less light curing contraction. **MASTER FLOW** presents high fluidity, provides access into microcavities. It is easy to handle and is not soluble in the oral environment. It is compatible with all adhesives and resins. It has radiopaque characteristics. The light curing time is 20 seconds.

TECHNICAL INFORMATION:

MASTER FLOW is a light cured, hybrid composite resin with organically modified glass-ceramic microparticles, giving the product greater resistance, greater shine, and less light curing contraction. **MASTER FLOW** presents high fluidity, provides access into microcavities. It is easy to handle and is not soluble in the oral environment. It is compatible with all adhesives and resins. It has radiopaque characteristics. The light curing time is 20 seconds.

USE INSTRUCTIONS:

- I - Definitive esthetic restorations as a initial layer:
1. Make the crown polishing with pumice powder and water;
 2. Accomplish the clinical procedures of cavity preparation according to required technique;
 3. Select the color with hydrated teeth;
 4. Place rubber dam;
 5. Apply the 37% Orthophosphoric Acid (a.e. Ataque Gel) during 15 seconds for dentin and 30 seconds for enamel.
 6. Rinse thoroughly with water for 30 seconds;
 7. . Dry the dentin (without dehydrate) with slight intermittent air jets or by capillarity. It is important to dry the enamel well and not to dehydrate the dentin;
 8. Apply the bonding agent (a.e. Master Bond DE) according to manufacture’s instructions;
 9. Apply MASTER FLOW to fill the cavity as the initial layer (thickness min. 1 mm and max 2 mm);
 10. Light Cure;
 11. Place with composite (a.e. Master Fill), according to manufacturer’s instructions.
- II - Esthetics dental restorations:
1. Make the coronal polishing with pumice powder and water;
 2. Accomplish the clinical procedures of cavity preparation according to required technique;
 3. Select the shade with the tooth still moist;
 4. Place rubber dam;
 5. Apply the 37% Orthophosphoric Acid (a.e. Ataque Gel) during 15 seconds for dentin and 30 seconds for enamel;
 6. Rinse thoroughly with water for 30 seconds;
 7. Dry but do not dehydrate the dentin though soft airflow;
- It is important to dry enamel but do not dehydrate Dentin;
8. Apply the bonding agent (a.e. Master Bond DE) according to manufacture’s instruction;
 9. Apply MASTER FLOW to fill the cavity in thin layers (thickness min. 1 mm and max 2 mm);
 10. Light cure for 20 seconds;
 11. Because of its viscosity (high glaze) final polishing is not necessary unless the occlusal or enamel adjustments are necessary.
- III - Prosthesis repair and indirect restorations:
1. Carry out, if necessary, the masking of the metal by means of an opacifier (eg. Biocal);
 2. Apply the Porcelain Conditioner for 4 minutes, rinse and dry;
 3. Apply Silane all over the conditioned area and wait for one minute; apply new layer;
 4. Apply bonding system (a.e. Master Bond DE) according to the manufacture’s instruction;
 5. Place MASTER FLOW fulfilling cavity in thin layers (thickness of min 1 max 2 mm);
 6. Finish and polish.

SPECIAL CARES:

- Dentin pulp complex must be protect with an appropriate base before etching;
 - Eugenol based liners affect negatively in polymerization;
 - Discoloration can occur when cationic mouth rinse is used as well as bacterial plaque disclosing and chlorhexidine;
 - To obtain a satisfactory light curing, high quality equipments must be used. The distance to restoration surface should be inferior to 5 mm when light curing; follow rigorously the instructions for thickness of restorative and light exposure time. Composites do not duly light cured can result a dentin pulp complex irritation, reducing the mechanical properties and discoloration.
 - When working on extensive cases, change the position of the operating light to avoid the early composite set-up.
- Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

STORAGE:

After using keep the syringes tightly closed to avoid exposure to the light. Do not store the product near a


product containing Eugenol. Eugenol interferences in polymerization.

CONTRAINDICATION:

MASTER FLOW should not be used on the impossibility to realize the absolute isolation. Should not be used in patients is known to have sensibility to any of the ingredients listed.

WARNINGS:

Do not use material after expiry date.
The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.
“KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.
FOR USE IN DENTISTRY ONLY

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