ideal for carrying and spreading the monomers over the surface. With their hydrophilic end, they are capable of efficiently moistening the dentin surface, while the hydrophobic part offers sufficient bonding points for copolymerization with the subsequent dentin adhesive. Dentin conditioners need to be removed after application. (João Galan Júnior, 1999).

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Perform absolute isolation.

I. Composite Restorations:

1. Etch the enamel for 30 seconds and the dentin for 15 seconds.
2. Rinse the conditioned area thoroughly with water/air for 30 seconds.
3. Dry the enamel completely; the dentin must be dried with an absorbent paper so that this tissue does not dehydrate, which could impair adhesion; to receive the adhesive application, the dentin must have a moisture shine.
4. Actively apply **MASTER BOND DE** evenly over the entire conditioned surface.
5. Evaporate the solvent with a jet of air.
6. Apply a new layer of **MASTER BOND DE** in the same way as in the previous step.
7. Evaporate the solvent with a jet of air;
8. Light cure for 20 seconds.
9. Apply the restorative material according to the manufacturer's instructions.
- II. In Indirect Ceramic, Metal and Composite restorations (Inlays, Onlays and Full Crowns):
1. Prepare the internal surface of the part, according to the material from which it was made.
2. Carry out prophylaxis on the tooth prepared with pumice stone and water.
3. Carry out acid etching (e.g. ATTAQUE GEL) for 30 seconds on enamel and 15 seconds on dentin.
4. Rinse the conditioned area thoroughly with water/air for 30 seconds.
5. Dry the dentin (without dehydrating it) with absorbent paper. It is important to dry the enamel well and not dehydrate the dentin; the dentin should have a moisture sheen.
6. Actively apply **MASTER BOND DE** evenly over the entire conditioned surface;
7. Evaporate the solvent with a jet of air.
8. Apply a new layer of **MASTER BOND DE** in the same way as in the previous step.
9. Evaporate the solvent with a jet of air;
10. Light cure for 20 seconds.
11. Use resin cement according to the manufacturer's instructions.

III. Repair of Metal Ceramic Restorations:

1. Roughen the exposed surface of the metal and ceramic with a jet of aluminum oxide or diamond burs.
2. Bevel the margins of the fractured porcelain.
3. Wash and dry.
4. Opacity the metal (e.g. Biocal).
5. To increase adhesion, apply porcelain conditioner for 4 minutes, wash and dry.
6. Apply Silane to the fractured porcelain and wait 1 minute; carry out a new application.
7. Apply **MASTER BOND DE** evenly over the entire conditioned surface.
8. Evaporate the solvent with a jet of air.
9. Apply a new layer of **MASTER BOND DE** in the same way as in the previous step;
10. Evaporate the solvent with a jet of air;
11. Light cure for 20 seconds.
12. Complete the repair with composite following the manufacturer's instructions.

*Mask the metal if necessary.

IV. Cementation of Porcelain Veneers:

1. Prepare the internal surface of the part, according to the material from which it was made.
2. Apply Silane to internal surfaces and wait for 1 minute; carry out a new application.
3. Perform preparation prophylaxis with pumice stone and water.
4. Wash and dry.
5. Carry out acid etching (e.g. ATTAQUE GEL) for 30 seconds on enamel and 15 seconds on dentin.
6. Rinse the conditioned area thoroughly with water/air for 30 seconds.
7. Dry the dentin (without dehydrating it) with absorbent paper. It is important to dry the enamel well and not dehydrate the dentin; the dentin should have a moisture sheen.
8. Actively apply **MASTER BOND DE** evenly over the entire conditioned surface.
9. Evaporate the solvent with a jet of air.
10. Apply a new layer of **MASTER BOND DE** in the same way as in the previous step.
11. Evaporate the solvent with a jet of air;
12. Light cure for 20 seconds.
13. Spatulate the cement and apply it to the veneers.
14. Carefully place the veneers and remove all excess.
15. Light cure first on the gingival margin followed by the crown and incisal.

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MASTER BOND DE is contraindicated for people who are sensitive to any of the product's components.

PRECAUTIONS:


- Handle the material with latex gloves to avoid contact with the skin.
- If contact occurs with the eye area, wash thoroughly with running water. If irritation persists, consult a doctor.
- The risk of postoperative sensitivity can be minimized by performing conservative cavity preparation.

- The restorative material must be placed and light-cured in layers to reduce polymerization contraction and occlusion must always be checked.
- When using the light curing device, follow the manufacturer's guidelines regarding eye protection.
- After use, the bottles must be perfectly sealed to prevent the material from being exposed to light.
- Packaging opened for a long time can lead to changes in the characteristics and properties of the product.
- Do not store near products containing Eugenol, as this interferes with polymerization.
- Dispose of the product in accordance with local legislation, ensuring its complete decharacterization in order to prevent its reuse and environmental damage.

WARNINGS: Do not use products beyond the indicated validity period.

The manufacturer is not responsible for damages caused by incorrect use or use not foreseen in the instructions for use.

KEEP THE PRODUCT AWAY FROM CHILDREN.
EXCLUSIVE USE BY A DENTAL SURGEON.

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MASTER BOND DE Light-cured Single-component Adhesive for Enamel and Dentin

PRESENTATION:
Introductory Kit:
1x5mL o 2x5mL bottle of MASTER BOND DE
1x5mL ATTAQUE GEL bottle;
20 microbrushes;
Replacement:
1x4mL and 1x5mL bottle of MASTER BOND DE
ANVISA nº 10298550091

COMPOSITION:

MASTER BOND DE: Bisphenol A Glycidyl Methacrylate, Dimethacrylate Groups, Alcohol, Sodium Fluoride and Catalysts.

INDICATION:

MASTER BOND DE is used in light-cured composite resin restorations to provide adhesion of the restoration to the tooth structure (dentin or enamel); in dentin desensitization and fragment bonding processes; as a complement to adhesive cementation systems; for bonding porcelain veneers and porcelain repairs.

TECHNICAL INFORMATION:

MASTER BOND DE is a dentin-enamel adhesive, which allows the bonding of the dental structure (enamel and dentin) to resino- materials. It is an adhesive system in which the primer and adhesive are found in a single bottle (single component), which facilitates the handling by the dentist. It is light-cured and easy to apply, it is hydrophilic, provides adhesive strength and protection against micro-leakage. An adhesive system, such as **MASTER BOND DE**, micro-mechanically joins the dentin to the resinous material by forming a layer of intertubular dentin impregnated by resin, known as a hybrid layer. Master Bond DE has low viscosity, which allows it to penetrate into the porosities of the dental tissue and mechanically imbricate it to obtain high adhesion. However, for adequate adhesion, it is necessary to chemically condition the enamel and dentin, with a 37% phosphoric acid, to completely remove the smear layer (smear layer that is deposited on the dentin tissue after cavity preparation), previously to the application of the adhesive. Primers are substances that contain hydrophilic and hydrophobic bifunctional monomers and are dissolved in alcohol which, as it has the property of displacing water and volatile properties, is