II - Metal Opacifier in Prosthesis Repair:

1. Roughen the exposed metal and ceramic surface with aluminum oxide jet or diamond burs;
2. Bevel the fractured porcelain margins;
3. Rinse and dry;
4. Select **BIOCAL** shade (dentin or white) and apply a thin layer (max. 0,2mm) directly into the cavity with the applicator nozzle, exclusively on the metal;
5. Light cure during 20 seconds;
6. Apply over **BIOCAL** the restoring material following manufacturer's instructions.

SPECIAL CARES:

After use keep syringes tightly close to avoid exposure to light.

Do not store **BIOCAL** near products containing Eugenol, because it interferes in the polymerization.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:

BIOCAL should not be used in patients is known to have sensibility to any of the ingredients listed.


WARNINGS:

Do not use material after expiry date.

The manufacture is not responsible for any damage resulting from a misuse or not provided for the instructions of use.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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BIOCAL
Light Cured Calcium Hydroxide
biodinâmica®

CONTENTS:

Reposition:
1x2g syringe BIOCAL Dentin.
3 application tips.
or
1x2g syringe BIOCAL White.
3 application tips.
Reg. ANVISA: 10298550061.

COMPOSITION:

Calcium, dimethacrylate groups, inorganic filler, photoactivator, stabilizer and pigments.

INDICATION:

BIOCAL is indicated for pulp protection and covering under restorative materials, cements and other base material and as a metal and dentin opacifier.

TECHNICAL INFORMATION:

BIOCAL is a photopolymerizable and radiopaque Calcium Hydroxide composition. **BIOCAL** no inhibits the polymerization of self- or light-cured restorative resins. In the dentine version, has a color similar to dentin, does not interfere with the final color of restorations in composite resins.

BIOCAL was formulated to present high resistance to dissolution by the phosphoric acid used in the acid etching of dental structures. Its formulation also offers potential adhesion to metal surfaces, acting as an opacifier in prosthesis repairs.

USE INSTRUCTIONS:

I - Indirect Pulp Covering and Reactional Dentin Opacifier:

1. Wash and dry the cavity;
2. Select **BIOCAL** shade (dentin or white) and apply a thin layer (max. 0,2mm) directly into the cavity with the applicator nozzle;
3. Remove material excess;
4. Light cure during 20 seconds.