
MASTER CEMENT DUAL

Cimento Resinoso

biodinâmica®

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6. Provide the etching of the dental structure with phosphoric acid 37% (ATTAQUE GEL) and leave it for 15 seconds on the dentin and for 30 seconds on the enamel. Then rinse thoroughly and dry without dehydrate the dental surface (a light glow of humidity should be observed on the surface).

7. Apply conventional adhesive agent according to technical indicated by material in use (for example MASTER BOND DE)

8. APPLICATION:

Using Base and Catalyst in double syringe and mixing tips:

8.1. Adapt the mixing tip to the syringe, and apply a thin layer of cement inside the prosthetic part or over the dental cavity covering the entire surface.

Using Base and Catalyst in separate syringes:

8.2. Dispense the required amount of **MASTER CEMENT DUAL** in a mixing block in the proportion of 1: 1 of Base paste and Catalyst paste. Mix the product for approximately 30 seconds and bring the cement to the restoration using a brush or spatula.

*Observation: **MASTER CEMENT DUAL**: Working time about 1,5 to 2 minutes in ambient temperature; Setting time about 2 to 3 minutes. The increasing of temperature accelerates the setting time. As it is a duo polymerization material, avoid its manipulation long before use, because excess of luminosity may accelerates setting time.

9. Put the restoration in position keeping it with slight pressure and remove excesses aided by a brush or an adequate instrument.

10. Light cure by 40 seconds in each segment, starting the proximal margins.

11. Remove light cured material excesses using specifics tips and realize the finishing with silicon tips and polishing disk.

SPACIAL CARE:

Eliminate completely all cementation proof material (MASTER CEMENT TRY IN) before do the definitive cementation. Residues of this material can damage the adhesion between the restoration and the resin cement;

Protect the pulp with a lining material or cavity base, if necessary.

Do not use materials containing Eugenol, because Eugenol influences negatively at the polymerization time.

Observe that exists a slight difference between the colors of MASTER CEMENT TRY IN and the correspondent colors of MASTER CEMENT LC, before the photo activation, due to the fact that the coloration is changed after the light curing. The color of MASTER CEMENT TRY IN is correspondent to the color of MASTER CEMENT LC, when completely light cured.

Packaging opened for a long time can lead to changes in product characteristics and properties.

Prohibited reprocessing.

CONTRAINDICATION:

The products components of the adhesive cementation system of **MASTER CEMENT DUAL** are contraindicated for people with related sensibility to any of its components

WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

“KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.

FOR USE IN DENTISTRY ONLY

or ceromers, inlays, onlays, crowns and bridges with or without metal, endodontic core, pre-fabricated or welded and adhesives prosthetics.

COMPOSITION:

Base Paste – **MASTER CEMENT DUAL:** Bisphenol A Glycidylmethacrylate; Dimetacrylates groups; Inorganic files; Catalysts; Stabilizers; Pigments.

Catalyst Paste – **MASTER CEMENT DUAL:** Bisphenol A Glycidylmethacrylate; Dimetacrylates groups; Inorganic files; Catalysts; Stabilizers; Pigments.

TECHNICAL INFORMATION:

Resin cements have the same basic composition as a composite resin, that is, a matrix of Bis-GMA (Bisphenol-A-glycidylmethacrylate) in combination with low viscosity monomers, plus inorganic material load, as well as curing initiators and pigments . They should be used after applying 37% phosphoric acid etching to the dental structure for 15 to 30 seconds and dentin adhesive. Resin cements can be divided into three groups according to the activation method: self cured, light cured and double-cured or dual resin cements, which involve chemical and light activation.

MASTER CEMENT DUAL is a dual adhesive cementation system, or double cured, which involves chemical and light activation, and allows the professional to choose the most appropriate cementing agent for each clinical procedure, guaranteeing the curing of the product in the presence or absence of light, in areas where exposure to light is critical. The product is indicated for adhesively bonding indirect restorations to the tooth structure.

It is compatible with conventional dentin/enamel bonding agents on the market, such as multiple-use adhesives as well as single-component adhesives (e.g. Master Bond DE). Its ceramic load gives the product excellent physical qualities such as high flexural strength and the combination of monomers provides a high degree of conversion during light curing. It is very easy to apply due to the presentation of the two components in a double-body syringe, ensuring application in appropriate proportions and homogeneity of the pastes, preventing the incorporation of bubbles in the product. Comes with mixing tips for applying the product.

Presented in different colors, which provides excellent aesthetic properties in the cementation of metal-free prostheses. In light curing, the indicated light curing time is 40 ± 10 seconds, using a light curing device with a minimum power of 400MW, with no restrictions on the light curing device used.

MASTER CEMENT TRY-IN, which may follow up the product, it is a try-in material for the cementation, it has glycerin base and is water-soluble. It acts promoting a simulation of the color that will be obtained through using

MASTER CEMENT DUAL and **MASTER CEMENT LC** (resin cement light curing), when completely polymerized. It is introduced in the same colors as resin cement **MASTER CEMENT DUAL** and **MASTER CEMENT LC**.

MASTER CEMENT LC, which may follow up the product, is a system of adhesive cementation light curing. This product is indicated to join indirect restorations till 2,0 mm to dental structures. It is compatible with conventional adhesives agents to dentin/enamel existents at market, as adhesives multiple use (for example Master Bond Multi-Use) as well with adhesives singles components (for example Master Bond). Your ceramic power gives to product greats physics qualities as elevated flexural resistance and the monomer combination gives high conversion grade in the light curing. The time for light curing indicated is about 40 to 10 seconds, using an equipment curing light with 400MW of minimum power, there is no restrictions about the equipment curing light used. It is introduced in many colors, what gives to it an excellent esthetics properties at cementation of facet and prosthesis without metal.

The acid present in the ATTAQUE F - (Porcelain/Ceramic Conditioner), which may follow up the product **MASTER CEMENT DUAL**, promotes a pattern of deep penetration and extensive of surface, making it micro-retentive, promoting a greater resistance of interface union composite porcelain/resin, resin porcelain/cement.

The Silane which may be in the **MASTER CEMENT DUAL KIT**, acts as adhesivity promoter between ceramic and composite, through a stable chemistry bond, providing an resistance increasing of ceramic/composite interface. It has chemical adhesion to ceramic silex, to resin cement and to matrix of composites resins..

INSTRUCTIONS OF USE:

1. Remove the temporary element and clean the tooth structure with pumice and water. In the case of root canals, clean with E.D.T.A. Trisodium Gel 24%, for removing the “sealer layer”, facilitating the penetration of the adhesive agent (e.g. Master Bond DE) into the dentinal canalculi, providing greater cementation resistance.

NOTE: Only use eugenol-free temporary cements (e.g. PROVITEMP) as eugenol can inhibit the curing of **MASTER CEMENT DUAL**.

2. Promote the isolation of the area to be worked.

3. Prove the restoration with relief of try-in paste (**MASTER CEMENT TRY-IN**). Evaluate which color is the most indicate for cementation.

4. Provide glycerin gel cleaning with water and dry with air jet water and oil-free.

5. Realize the treatment of the restoration according to indicate technical by the material in use. Porcelain restorations and composites must be treat in them internal face with ATTAQUE F (porcelain conditioner) for 1 minute. Then wash and dry with air jet water and oil-free, and realize the SILANE application. Wait for 1 minute and repeat the application.

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Resin Cement

biodinâmica®

1x5g Double Syringe **MASTER CEMENT DUAL** (Auto and Light Curing); + accessories for exclusive use with the product; or

1x2,5g syringe base paste **MASTER CEMENT DUAL** (Auto and Light Curing) proportion 1:1; and 1,2,5g syringe catalyst paste **MASTER CEMENT DUAL** (Auto and Light Curing) proportion 1:1+ accessories for exclusive use with the product.

Available in the colors: TRANSLUCENT - A1 - A2 - A3 - W.B. (White Bleaching) - W.O. (White Opaque) White Opaque.

ANVISA nº 10298550114.

INDICATION:

MASTER CEMENT DUAL is a resin cement, indicated to adhesive cementation of indirect ceramic restorations, composite