
# HIDRÓXIDO DE CÁLCIO P.A

## Capeedor Pulpar

biodinâmica®

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1x10g flask **CALCIUM HYDROXIDE P.A.**  
ANVISA n° 10298550086.

### INDICATION:

Intracanal curative material. Excellent use on stimulating the repairing dentin formation (increasing its thickness and giving more protection to the pulp tissue).

### COMPOSITION:

Calcium Hydroxide P.A. Powder (99 to 100,5%).

### TECHNICAL INFORMATION:

**CALCIUM HYDROXIDE P.A.** is now one of the most accepted material to induce the formation of the reparative dentine, presenting indispensable properties for materials that act in procedures on dentine. These properties are:

a) It allows an increase of the calcification of the layer of pre-existing dentine giving better protection conditions to the pulp when submitted to lesions;

b) It avoids that the procedure or the used material causes irritation to the adjacent structures;

It is used in Endodontic as intracanal material acting as a potent antimicrobial by the inhibition of the bacterial enzymes;

It has excellent histocompatibility and it stimulates the host tissue, propitiating the tissular repair (for tissular enzymatic activation, creating mineralization effect).

### USE INSTRUCTIONS:

I - To prepare a 0,2% Calcium Hydroxide Solution:

1) In an appropriate recipient, add 100 ml of distilled water to the content of the flask of the product (10g of **CALCIUM HYDROXIDE P. A.**);

2. Maintain this mixture resting, so that the excess of **CALCIUM HYDROXIDE P.A.** stays deposited in the bottom of the recipient.

NOTE:

- It is not necessary to shake for the use, because the alkaline solution is above the deposition.

- This alkaline solution should be used before pulp protection and restoration are placed, in all cavity types, any deep.

- **CALCIUM HYDROXIDE P.A.** deposited in the bottom of the recipient can be used by professional, because it constitutes an excellent paste for direct protections.

II - In Endodontic:

1. Mix **CALCIUM HYDROXIDE P.A.** with the vehicle of your choice (distilled water, Camphorated

Paramonochlorophenol or others);

2. Fill the canal with the obtained paste.

### SPECIAL CARES:

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

### CONTRAINDICATION:

**CALCIUM HYDROXIDE P.A.** should not be used in patients is known to have sensibility to the product.

### WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.  
"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

FOR USE IN DENTISTRY ONLY