BIODAM should not be used in patients known to be sensitive to any of the listed ingredients.

It should not be used in the following cases:

- Patients who already have gingivitis or periodontal illness and gingival injuries;
- Patients who have undergone recent periodontal surgery;

WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not in accordance with these instructions.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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
Reg. Anvisa nº 10298550101

COMPOSITION:

Dimethacrylate groups; silicon dioxide; BHT; catalysts; inorganic filler and pigments.

INDICATION:

BIODAM is indicated as a gingival protector in cases of use products that pose risks of irritation are used (e.g.: bleaching agents based on hydrogen peroxide for clinical use), preventing the contact of the irritating agent with the gingival tissue.

TECHNICAL INFORMATION:

For the professional bleaching procedures with hydrogen peroxide or for the gingival isolation from any other product, is necessary to use a material that offers protection. **BIODAM** is an insulating and protective material. After being light-cured creates a rigid mechanical barrier to isolate the teeth and protect the gingival tissues. BIODAM can substitute for or complement the rubber dam for the gingival protection.

BIODAM promotes a relative isolation, adhering to the gingiva and mechanically blocking the passage of any product crossing. It has easy application. The colors facilitate the contrast with the tooth. It adheres well to the gingiva, and its removal is easily done with the help of an instrument (e.g.: explorers), detaching the product.

USE INSTRUCTIONS:

1. Thoroughly dry the teeth and the gingiva, especially the gingival portion;
2. Apply a thin layer of **BIODAM**, approximately from 4 to 6 mm of height and from 1 to 1,5 mm of thickness on the gingiva. Cover about 0,5 mm of the enamel, in the cervical region, to prevent the contact with the gingival border and repeat the process in all teeth that will be isolated.

Comment: In the presence of diastemas, also expand the lingual layer, and protect the mesiodistal.

Extending **BIODAM** between the teeth provides greater retention of protection.

3. Light cure for 20 to 30 seconds and continue this procedure until the entire length of **BIODAM** has been light cured.

4. Verify if the product is adequately light-cured and if it offers protection on the entire gingivas surface needed.

5. For removal, detach the product carefully with the help of an instrument.

SPECIAL CARES:

Light cure the product in accordance with the use instructions, respecting the time and the irradiation site.

Incomplete or excessively long curing may cause slight gingival irritation.

Dispose of the product in accordance with local regulations, ensuring its complete disposal, preventing its reuse and minimizing environment damage.

CONTRAINDICATION: