7. Protect the patient, assistant and professional with safety equipment (gloves, goggles, cap, mask and apron) to handle the product;
8. To prepare the gel, uncover the syringes, one from the HYDROGEN PEROXIDE phase and the other from the THICKENING phase, and connect them directly to each other. The syringes have been specially designed to be connected without the need for a connector. Make sure the syringes are securely connected. Make the preparation away from the patient;
9. Push the plungers back and forth 10 times or until the product reaches a homogeneous color, indicating complete mixing and activation of the gel.
10. After activation, the gel must remain within the HYDROGEN PEROXIDE phase syringe so that the applicator tip can be connected. Disconnect the syringes and connect the applicator tip;
- NOTE: The applicator tip does not fit into the THICKENING phase syringe.
11. Apply a thin layer of the gel to all the teeth to be treated. The thin layer will prevent the gel from flowing.
- Caution: Do not lift the tip when changing teeth. Make continuous movements. Make sure all dental surfaces are completely covered by the gel;
12. The gel should remain in contact with the surface to be clarified for 45 consecutive minutes without the need for reapplication in the same session. In patients with a history of output sensitivity, one can reduce the time to 30 minutes of stay and increase the number of visits;
13. At the end of the bleaching procedure, suction the gel excess, remove with gauze and wash with plenty of water, being careful to suck all the water to avoid the intake of the product;
14. If desired, the GHF desensitizer may be applied without damage to bleaching. The application should be made before removal of the gingival barrier.
15. Remove the gingival barrier by detaching it forward;
16. If you do not apply the GHF, it is recommended to apply a fluorine-containing material (e.g. FRUTTI-FLUOR GEL NEUTRO) for 1 to 4 minutes for the mineralization of demineralized areas;
17. Compare with the color initially selected.

Obs.: If there is gel remaining in the syringe, do not use it on another patient, discard it because after mixing the HYDROGEN PEROXIDE and THICKENING phases, the gel should be used immediately to maintain the whitening action.

SPECIAL CARE:

- BEFORE BLEACHING:

1. Assess and determine patient expectations before starting treatment. Inform the patient that pre-existing restorations, crowns, and prosthetics do not undergo bleaching agent action. Inform that due to unforeseen factors or tooth enamel variations, some teeth may not meet expectations;
2. Examine the patient and ensure proper oral health conditions such as: healthy gingival tissues, plaque control, caries-free teeth, infiltration of restorations, or other conditions that may compromise bleaching treatment.
- 2.1. Carry out the removal of calculus, extrinsic stains and other superficial pigmentations before bleaching consultation to avoid any bleeding that may interfere with the application and adaptation of the protection barrier and whitening treatment.
- 2.2. If you use sodium bicarbonate for prophylaxis, start the bleaching treatment only after one week.
3. Treating dental hypersensitivity, if present, as it is common to exacerbate the painful symptomatology during bleaching treatment.
4. Register the color of the teeth by means of color scale and/or photographs, evaluate and inform the patient clearly about the results that can be expected. Start from the top arch.
5. Guidance on possible inter-relationships:
- 5.1. Teeth sensitivity to cold, which disappears 1 to 3 days after discontinuation or termination of treatment;
- 5.2. Teeth that have restorations using large amounts of amalgam may appear darkened after the bleaching treatment. This occurs because the amalgam appears more visible through the whitened enamel. For these cases, the replacement of amalgam restorations with composite resin (e.g. MASTER FILL) after tooth whitening is recommended.
6. Do not anesthetize teeth for whitening so as not to mask patient sensitivity.
7. Clinically INVISIBLE hypocalcified areas prior to bleaching may become visible afterwards. Because of the differences in mineralization between these areas and the normal enamel, they brighten before and momentarily become visible. Continue with bleaching treatment and reevaluate for 7 to 14 days. In the vast majority of cases, the continuity of bleaching blends this color difference, homogenizing it, after total rehydration of the teeth. The use of fluoride may help the disappearance of these areas. If there is a VISIBLE color difference before the beginning of the treatment, the lighter areas can be covered with BIODAM so that the final result is homogeneous.
8. Pre-existing restorations and crowns do not undergo bleaching agent action and cannot be cleared. Thus, whitening treatment should be performed prior to replacement of aesthetic restorations.

- AFTER BLEACHING:


1. Check the presence of sensitivity and, if appropriate, treat appropriately using the desensitizer indicated as the case may be. (e.g. FRUTTI-FLUOR GEL, G.H.F.).
2. Encourage the patient to maintain the bleaching effect obtained by the treatment, returning to the revision visits, for evaluations and strengthening sessions if necessary.
3. If it is necessary to replace existing aesthetic restorations, it is recommended to wait 7 to 15 days after bleaching, as during this period the color stabilization occurs and the elimination of any traces of oxygen derivatives, which may interfere with the polymerization adhesive systems, composite resins and resin cements used in restorative treatment.

CONTRAINDICATION:

- Already installed tables of gingivitis, periodontal disease and gingival lesions;
- Patients who underwent periodontal surgical intervention recently;
- Pregnant women and women who are breastfeeding;
- Patients with cancer;
- Patients allergic or sensitive to any of the components of the product.

WARNINGS:

Exclusive use of dentist.
Do not use material after expiry date.
The manufacture is not responsible for any damage resulting from a misuse or not provided for the instructions of use.
SELLING UNDER PRESCRIPTION
"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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CLARIDEX OFFICE 35%

Dental Whitening Agent (Hydrogen Peroxide)

biodinâmica

CONTENTS:

- 1 syringe set of 1,8 g being 1 thickener syringe with 0,42 g and 1 hydrogen peroxide syringe with 1,38 g + accessories.
- or
- 3 syringe set of 1,8 g being 3 syringe of thickener with 0,42 g and 3 syringe of hydrogen peroxide with 1,38 g + 1 x 2 g of syringe of Biodam + accessories.
- or
- 6 syringe set of 1,8 g being 6 thickener syringe with 0,42 g and 6 hydrogen peroxide syringe with 1,38 g + 1 x 2 g Biodam syringe + accessories.
- ANVISA nº 10298550100

INDICATION:

CLARIDEX OFFICE 35% is indicated for professional dental whitening, for exclusive use in dentistry.

COMPOSITION:

Hydrogen peroxide Phase: hydrogen peroxide and thickener.
Thickening Phase: thickener, preservative and activator.

TECHNICAL INFORMATION:

CLARIDEX OFFICE 35% is a dental bleaching based on hydrogen peroxide, composed of 2 phases: hydrogen peroxide phase and thickener phase. When the phases are mixed together they form a gel, which results in the concentration of 35% hydrogen peroxide. The pH, after mixing the phases, is between 6.0 and 7.0 and remains stable until the end of the procedure, which avoids demineralization of the dental surface and does not generate sensitivity. Bleaching agents act on the tissues by oxidation and/or reduction of the darkening elements, their action is related to the ability to penetrate the dental structures (enamel and dentin). There is a reaction between the agent and the pigmented molecules, which are oxidized and become simpler and less pigmented, causing the pigments incorporated into the tissues to be totally or partially eliminated from the dental structure by diffusion. Most stains are susceptible to chemical oxidation from the bleaching procedure, however, the older the staining, the more difficult it can be to remove them. Dentin staining caused by the presence of tetracycline salts in general, oxidation requires more time for action and the result will depend on the degree of severity. Spots from metallic pigmentation (amalgam and certain endodontic sealing cements) are difficult to remove.

CLARIDEX OFFICE 35% presents ideal viscosity for a better application and should be used exclusively by a professional in a dental practice, as it requires the use of an efficient intraoral tissue protection system (e.g. gingival barrier, e.g. BIODAM or absolute insulation), since the peroxide oxidizing action can damage the soft tissues.

It can be used alone or in conjunction with the self-applied bleaching system. It is not necessary to apply a light source as an additional method to accelerate the bleaching process; the product will act effectively in the absence of light curing because its action is independent of the presence of light. It is suggested that the product remains in contact with the dental surfaces for a period of 45 consecutive minutes, there is no need for sectioned application every 15 minutes. This fact is due to the properties of the product keep its neutral pH and action continues for a longer period of time after mixing. For best results other applications may be required in subsequent consultations, depending on the case and response to bleaching treatment. In patients with a history or sensitivity potential, one can reduce the time to 30 minutes of stay and increase the number of visits. In addition, the GHF desensitizing agent can be applied immediately after the removal of **CLARIDEX OFFICE 35%** without damage to the bleaching process.

INSTRUCTIONS FOR USE:

1. Assess the clinical condition of the patient in order to detect if any previous treatment is necessary;
2. Perform a prophylaxis with pumice and water (do not use sodium bicarbonate);
3. Position the labial retractor;
4. Register the color of the patient's teeth by means of a color scale before starting bleaching;
5. Perform relative isolation using a light-curing gingival barrier, e.g. BIODAM. Apply a thin layer of BIODAM approximately 4 to 6 mm in height approximately 1 mm thick over the gingiva. Cover about 0.5 mm of the enamel, in the cervical region, to avoid contact with gengivla. Repeat the process on all teeth to be bleached;
6. Make sure that all gums around the teeth have been covered and if the product is fully polymerized and protecting the entire surface of the gum. Note: In the presence of diastema, extend the layer to the lingual, and protect the niches;