of the material.

SPECIAL CARES:
Both, Liquid and Cement, when in contact with the skin and mucosa are corrosive.
KEEP LIQUID DISTANT FROM THE EYES AREA AND AVOID LONG CONTACT WITH SKIN AND ORAL TISSUE.
In case of contact with oral tissue and skin, remove it immediately with cotton embbeded in alcohol and wash with water.
In case of contact with the eyes, wash immediately with water and seek medical advice.

SPECIAL RECOMMENDATIONS:
Shake the powder before using;
Liquid flask must be carefully manipulated avoid any air gaps;
Use a rigid manipulating spatula and the plate area should be as small as possible;
The material surface must be always shining after spatulation;
Insertion can be done using a spatula, application tips or insertion syringes;
Condensation is made with condensers, burnishers or small ball of polyethylene foam.
Strengthen the material against the walls for bubbles removal and better adaptation.
Superficial protection must be realized owing to property that has this material to lose water at the beginning and incorporate after some time, that would cause the loss of mechanic propriety.
Right after the use of the manipulation instruments, rinse them with cold water to avoid solidification and a difficult removal of the residual material;
Do not use liquid and powder from different ionomers flasks;
After finish any of **BIOGLASS R** components both flaks should be disposed;
Avoid direct contact with EUGENOL based products once it interferes in some glass ionomer's working time.
The use of a cold plate increases material working time.
Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

STORAGE:
Do not keep in refrigerator;
Keep the liquid and powder flasks tightly closed;
After use close the products immediately;
Do not expose product to intensive light;
In case the liquid becomes more viscous, immerse it on water at 75°C temperature during 10 to 15 minutes to be reused.

CONTRAINDICATION:
BIOGLASS R should not be used in patients is known to have sensibility to any of the ingredients listed.


WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

“KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.

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COMPOSITION:

BIOGLASS R POWDER: Calcium, Barium and Aluminum Fluorsilicate, Lyophilized Polyacrylic Acid, Barium Sulfate and pigments.

BIOGLASS R LIQUID: Tartaric Acid and Deionized Water.

INDICATION:

BIOGLASS R is indicated for restorations class I,III and V; restorations class I and II in deciduos teeth; and in restorations of erosion, abrasion and radicular decay;
To rebuild cores where there is a still standing coronary structure (50%);
Fulfilling of small retentions;
Specially indicated for ART technique (atraumatic restorative treatment) in odontopediatrics, special patients and geriatric;
In the lining of cavities including light cured resins (sandwich technique);
For fissure sealing;
For recovery of cusped rigidity without dentin support.

TECHNICAL INFORMATION:

The glass ionomer is a prominent material in modern dentistry, mainly because its adhesion property and fluorine release. Such properties guarantee an excellent cavity sealing and decay prevention. Easy to manipulate. Excellent aesthetics and mechanical properties. Used in reconstruction of prosthetic pieces only in elements with more than 50% of dental structure. Fluor release guarantee an anticariogenic action. Continuous recharged of Fluor considering its bucal absorpition. It allows both the remineralization of adjacent tissues and the recovery of rigidity of cusps without dentin support. It has chemical adhesion to the dental structure through the connection with calcium ions, Higher adhesion to the enamel than the dentin. Similar thermal expansion and contraction to dentin, having no post-operating sensitivity. High adhesivity material having no need of additional retentions. It is biocompatible and hydrophilic. It´s the most compatability material to our dental structure. In spite of initial cement acidity, Calcium Hydroxide based cement is recommended for deep cavities. Fluor contributes to reduce material pH. Radiopaque material and the high adhesivity to the dental structure avoid marginal leakage.

USE INSTRUCTIONS:

PREPARATION OF CAVITY:

1. Prepare the cavity;
2. Clean cavity with water and dry, but not dehydrate;
3. In case of pulp being exposed, recommended to apply Calcium Hydroxide cement (a.e. BIOCAL);
4. Using a brush apply a dentin conditioner (a.e. Liquid BIOGLASS R) during 15 seconds to remove the smear layer. Wash with water and dry, but not dehydrate.

MIXING, WORKING AND SETTING TIME: (23 C and 50% UR):

Mixing time: 35 to 40 seconds/ Working time: 1 min. 30sec. to 2 min. 30 sec.

Setting time: 5 to 6 minutes.

MIXING RATIO:

Power / liquid ratio: 1 (one) measure of power to 2 (two) drops of liquid/deal ratio: (m/m): 1,6:1.

SHAKE LIQUID AND POWDER FLASK BEFORE USE.

For a precise dosage of the liquid turn the flask vertically down. Release the air from inside the flask.

After use, close it tightly both flasks to avoid humidity contact.

RESTORATIVE TECHNIQUE:

1. Remove remaining moisture from the cotton ball surface or carefully with the air jet.. Do not dehydrate;
2. Apply the cement into the cavity using a proper device. Not allow air gaps in the mixed material;
3. After adhesion of material, for 5 to 6 minutes apply a varnish or surface sealant .In case of excess, it must be removed after 15minutes. A varnish or surface sealant should be applied again. The final polishing must be realized 24 hours after the insertion