# Master Fill

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## Master Fill Light Cured composite for Restoration

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1x3g syringe Ataque Gel (37% Orthophosphoric Acid);

1x5mL flask Master Bond DE (Light Curing Single Component Adhesive for Enamel and Dentin);

Accessories.

Reposition:

1x4g syringe **MASTER FILL**

Available in shades: A1; A2; A3; A3,5; A4; B0,5; B1; B2; B3; B4; C1; C2; C3; C4; D2; D3; D4; Incisal; Incisal Blue; OA2; OA3; OB1; OB2; OC2.

Reg. ANVISA: 10298550040

### COMPOSITION:

Metacrylic monomer matrix; inorganic filler, initiators, stabilizers and pigments.

### INDICATION:

It is a micro hybrid resin indicated for anterior and posterior definitive dental restorations in enamel and dentin; splinting of mobile teeth; prosthesis repair; indirect restorations.

### TECHNICAL INFORMATION:

**MASTER FILL** is a high density composite resin for use in both anterior and posterior tooth restorations on enamel and dentin.

Its formulation contains inorganic particles varying between 0,04 and 2,2  $\mu\text{m}$ .

It presents a total inorganic filler volume of 79%.

**MASTER FILL** presents fluorescence properties, excellent finishing, polishing and mechanics properties, allowing its use at anterior restorations, veneers and occlusal surfaces. Setting time of 40  $\pm$  10 seconds, if using a 400 MW curing light (The indicated light curing range is 16 Joules).

### USE INSTRUCTIONS:

1 - Definitive esthetic restorations:

1. Perform prophylaxis with pumice and water;

2. Select the shade with the tooth still moist;

3. Do the cavity prepare in accordance with the chosen technique

4. Proceed the absolute isolation;

5. Do the etching (a.e. Ataque Gel), for 30 seconds at the enamel and for 15 seconds at the dentin;

6. Rinse thoroughly with water during 30 seconds;

7. Dry but do not dehydrate the dentin. It is important to dry well the enamel but do not dehydrate the dentin.

Observation: In case a self-etching adhesive system is used, in items 5 to 7, follow the manufacturer's

instructions.

8. Apply the bonding agent (a.e. Master Bond) as indicated by the manufacturer's instructions;

9. Apply **MASTER FILL** in small increments (maximum 2 mm) and light cure each one of the increments, until the desired anatomy is obtained;

10. Make the occlusal adjustment ;

11. Finish and polish in a subsequent session with sanding discs and abrasive silicone tips;

12. It is recommended the use of a surface sealant for restoration (a.e. Bioforty) according to manufacturer's instructions.

II - Prosthesis Repair and Indirect Restorations:

1. To improve the adhesion, use the porcelain conditioner according to the manufacturer's instructions;

2. Apply the Silane on the broken porcelain and wait for 1 minute, do a new application;

3. Apply bonding system (a.e. Master Bond) according to the manufacturer's instruction;

4. Apply **MASTER FILL** in increments of no more than 2 mm and light cure between increments;

5. Promote finishing and polishing in a subsequent session.  
Observation: Perform, if necessary, the masking of the metal.

### SPECIAL CARES:

- The dentin-pulp system must be protected with an appropriate base before the etching;

- Eugenol based liners affect negatively in the light curing;

- The use of cationic collutorium, plaque indicators and clorexidine may interfere at the coloration of the material.

- To obtain a satisfactory polymerization, adequate equipments are required and the distance for the polymerization of the surface must be inferior of 5 mm. Follow rigorously the light exposure time, as well as the thickness of the layer; Non polymerized composites may cause irritation to the dentin-pulp system, changes in the shades and also reduce the mechanical properties.

- Absolute isolation is recommended.

- The application of surface sealant is recommended (a.e. Bioforty), in order to improve the life time of the restoration.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

### STORAGE:

After using keep the syringes tightly closed to avoid exposure to the light.

Do not store near products containing Eugenol. Eugenol interferes in the polymerization.

### CONTRAINDICATION:


**MASTER FILL** should not be used in patients is know to have sensibility to any of the ingredients listed.

### WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

“KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.  
FOR USE IN DENTISTRY ONLY

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