RESIGLASS F

Light Cured Glass Ionomer Cement for LINING - A3

biodinâmica®

CONTENTS:

1 x 10g flask **RESIGLASS F POWDER**.

1 x 8mL flask **RESIGLASS F LIQUID**.

1 measuring spoon.

ANVISA nº 10298550066.

COMPOSITION:

RESIGLASS F Powder: Calcium, Barium and Aluminum Fluorosilicate, Polyacrylic Acid, Inorganic Fillers and Pigments.

RESIGLASS F Liquid: Dimetacrylates Groups and Catalyst.

INDICATION:

RESIGLASS F is a glass ionomer cement, hybrid (modified by resin), light cured, radiopaque, used for lining under restorations with resins, amalgam, ceramics and restorative metals.

TECHNICAL INFORMATION:

Glass ionomer is a very important material in the current and modern Odontology, mainly because its adhesion property and fluorine release. Such properties guarantee an excellent cavity sealing and decay prevention. Easy to manipulate. Excellent aesthetics and mechanical properties. Used in reconstruction of prosthetic pieces only in elements with more than 50% of dental structure. Fluor release guarantee an anticariogenic action. Continuous recharged of Fluor considering its bucal absorption. It allows both the remineralization of adjacent tissues and the recovery of rigidity of cusps without dentin support. Chemical adhesion to dental structures through its link with Calcium ions. Higher adhesion to the enamel than the dentin. Similar thermal expansion and contraction to dentin, having no post-operating sensitivity. High adhesivity material having no need of additional retentions. Biocompatible and hydrophilic material. It's the most compatibility material to our dental structure. In spite of initial cement acidity, Calcium Hydroxide based cement is recommended for deep cavities. Fluor contributes to reduce material pH. Radiopaque material and the high adhesivity to the dental structure avoid marginal leakage.

USE INSTRUCTIONS:

TEETH PREPARATION:

1. Prepare the cavity;
2. Clean cavity with water and dry, but not dehydrate;
3. In case of exposed pulp, recommended to apply Calcium Hydroxide cement (a.e. BIOCAL);
4. Pre treatment of dentin using **RESIGLASS F Liquid**, in the cavity, brushing it during 15 seconds. Wash with water and dry, but not dehydrate.

MIXING AND WORKING TIME: (23º C and 50% UR):

- Mixing time: 25 to 35 seconds

- Working time: 5 ± 0,5 minutes

MIXING RATIO:

- Power / liquid ratio: 1 (one) measure of power to 1 (one) drop of liquid

- Ideal ratio: (m/m): 1,8:1

SHAKE LIQUID AND POWDER FLASK BEFORE USE.

For a precise dosage of the liquid turn the flask vertically down. Release the air from inside the flask. After use, close it tightly both flasks to avoid humidity contact.

LINING TECHNIQUE:

- 1) Avoid contamination by water or saliva during application and insertion of the cement. Rubber dam is the best technique to obtain a perfect isolation.
- 2) Place spatulated cement on dental surface of cavity in a thin coat (0,5mm) using an adequate instrument;
- 3) Light cure during 20 seconds;
- 4) When a thicker coat is necessary, use the incremental coats technique, light curing each coat during 20

seconds;

5) To adhere to the composite, amalgam, ceramics or others, use a light cure Bonding material (a.e. MASTER BOND).

SPECIAL CARES:

- Both Liquid and Cement when in contact with the skin and mucosa are corrosive.
- KEEP LIQUID DISTANT FROM THE EYES AREA AND AVOID LONG CONTACT WITH SKIN AND ORAL TISSUE.
- In case of contact with oral tissue and skin, remove it immediately with cotton imbibed in alcohol.
- Rinse with water.
- In case of contact with the eyes, wash immediately with water and seek medical advice.

SPECIAL RECOMMENDATIONS:

- Shake the powder and liquid before using;
- Liquid flask must be carefully manipulated to avoid any air gaps;
- Use a rigid manipulating spatula, the plate area should be as small as possible;
- Material surface must be always shining after spatulation;
- Insertion can be done using a spatula, application tips or insertion syringes;
- Condensation is made with condensers, burnishers or styrofoam pellets;
- Reinforce material against the walls to remove any air gaps and obtain a better adaptation;
- Right after the use of manipulation instruments, rinse with cold water to avoid solidification and difficult removal of material residuals;
- Do not use liquid and powder from different ionomers flasks;
- Avoid direct contact with Eugenol based products once it interferes in some glass ionomer's working time.

After finish any of **RESIGLASS F** components both flasks should be disposed;

After use keep flasks tightly closed to avoid exposure to humidity and luminosity.

Avoid direct contact with Eugenol based products once it interferes in some glass ionomers' working time.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:


RESIGLASS F should not be used in patients is known to have sensibility to any of the ingredients listed.

Do not use material after expiry date.

The manufacture is not responsible for any damage resulting from a misuse or not provided for the instructions of use.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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