EPIC

Resina Composta Nano híbrida Fotopolimerizável

biodinâmica®

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1x4 g seringa de EPIC Corpo / Universal nas cores: A1; A2; A3; A3.5; A4; B1; B2; B3; B4; C1; C2; D2; D3; D4.
1x4 g seringa de EPIC Esmalte nas cores: EA1; EA2; EA3; EA3,5; EA4; EB1; EB2; EB3; EB4; EC1; EC2; EC3; EC4; ED2; ED3; ED4.

1x 4 g seringa de EPIC Dentina nas cores: DA1; DA2; DA3; DA3,5; DA4; DB0,5; DB1; DB2; DB3; DB4; DC1; DC2; DC3; DC4; DD2; DD3; DD4.

1x 4 g seringa de EPIC Dentes Clareados: WB; WO; WE; BL1 e BL2.

1x 4 g seringa de EPIC Opaquer na cor A1.

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EPIC

Light curing Nano hybrid Composite Resin

biodinâmica®

CONTENTS:

1x4 g syringe of EPIC Body / Universal in colors: A1; A2; A3; A3.5; A4; B1; B2; B3; B4; C1; C2; D2; D3; D4.

1x4 g syringe of EPIC Enamel in colors: EA1; EA2; EA3; EA3.5; EA4; EB1; EB2; EB3; EB4; EC1; EC2; EC3; EC4; ED2; ED3; ED4.

1x 4 g syringe of EPIC Dentin in colors: DA1; DA2; DA3; DA3.5; DA4; DB0.5; DB1; DB2; DB3; DB4; DC1; DC2; DC3; DC4; DD2; DD3; DD4.

1x 4 g syringe of EPIC Whitenend Teeth: WB; WO; WE; BL1 and BL2.

1x 4 g syringe of EPIC Opaquer in color A1.

1x 4 g syringe of EPIC Effects: Translucent and Incisal.

1x 2 g syringe of EPIC Pigments: Yellow; Blue; Brown; Ochre and Purple.

1x 2 g syringe of EPIC Flow in colors Opalescent and Opaquer B1.

ANVISA nº: 10298550210.

COMPOSITION:

Matrix of methacrylic monomers, inorganic fillers, initiators, stabilizer, preservative and pigments.

INDICATION:

EPIC is a nano hybrid resin indicated for definitive esthetic dental restorations anterior and posterior; fixing teeth with mobility, prosthesis repairs and indirect restorations.

TECHNICAL INFORMATION:

EPIC is a high-density, light curing, radiopaque and fluorescence nanohybrid composite resin. It contains in its inorganic composition a mixture of micro and nanoparticles, ranging in size from 5 nm to 2 µm. This size distribution allows for more efficient charge compaction, where the smaller particles fill the spaces between the larger particles. The inorganic composition is equivalent to approximately 79% by weight.

EPIC has fluorescence and excellent finishing and polishing, as well as excellent mechanical properties, which can be used both for previous restorations, including facets, and for posterior restorations. Complete color system, with colors for enamel, dentin, body, effects and pigments.

It has in its composition a monomeric organic matrix whose polymerization process occurs through the formation of free radicals. The camphorquinone present in the composition is activated by the incidence of visible light with a wavelength varying between 400 and 500 nm (blue spectrum) with maximum absorption at 468 nm. When activated, this component transfers one of its electrons to a supporting substance in the photoinitiation, resulting in the formation of free radicals and triggering the polymerization process.

Several factors contribute to an efficient light curing: incidence of light in the spectrum from 400 to 500 nm, power of the light curing unit, and distance from the light source, exposure time and thickness of the material layer.

The polymerization time for the EPIC resin is 40 ± 2 seconds if a curing light with a wavelength of 400 to 500 nm is used.

USE INSTRUCTIONS:

I - DEFINITIVE AESTHETIC DENTAL RESTORATIONS:

1. Select colors with the moistened teeth;
2. Promote coronary polishing with pumice and water;
3. Perform the cavity preparation according to the recommended technique;
4. Perform absolute isolation;
5. Perform the etching (e.g. ATTAQUE GEL) of the enamel for 30 seconds and dentin for 15 seconds;
6. Wash the etched area with plenty of water for 30 seconds;
7. Dry the dentin (without dehydrating) with light intermittent jets of air. It is important to dry the enamel well and not to dehydrate the dentin;
8. Use the adhesive system (e.g. MASTER BOND), following the manufacturer's guidelines;
9. Apply the resin in thin layers of 1 to 2 mm, light curing each layer;
10. Promote finishing and polishing;
11. Perform occlusion control;
12. Polish with sanding discs and abrasive silicone tips;
13. It is recommended to apply surface sealant for restorations (e.g. BIOFORTY), following the manufacturer's guidelines.

II - REPAIR OF PROSTHESES AND INDIRECT RESTORATIONS:

1. To increase adhesion, apply the porcelain conditioner, for the time recommended by the manufacturer, wash and dry;
 2. Apply SILANO over the fractured porcelain and wait for 1 minute; perform a new application;
 3. Use the adhesive system (e.g. MASTER BOND), following the manufacturer's guidelines;
 4. Apply the resin in thin layers of 1 to 2 mm;
 5. Promote finishing and polishing.
- * Perform, if necessary, the masking of the metal.

SPECIAL CARES:

- The dentin-pulp complex must be protected with an appropriate base before etching.
- Eugenol-based tools have a negative influence on the light curing time.
- The use of cationic mouthwashes, plaque detectors and chlorhexidine can interfere with the color of the material.
- For satisfactory light curing, use a light curing unit of adequate quality, and the distance for the light curing of the restoration surface must be less than 5 mm; carefully observe the exposure time, as well as the material thickness; Resins not properly cured can cause irritation to the dentin-pulp complex, as well as reducing mechanical properties and causing color change. In long-term work, deflect the light from the reflector to avoid premature hardening of the material.
- It is recommended to use absolute insulation.
- It is recommended to apply surface sealant for restorations (e.g. BIOFORTY), to seal the margins of restorations and increase the longevity of restorations.
- Dispose of the product in accordance with local legislation, ensuring its complete lack of characterization in order to prevent its reuse and environmental damage.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:

EPIC is contraindicated for people with reported sensitivity to any of the product's components.

WARNINGS:

- Do not use material after expiry date.
- The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.
- “KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.**
- FOR USE IN DENTISTRY ONLY**

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- 1 jeringa de 4 g de EPIC Cuerpo / Universal en los colores: A1; A2; A3; A3.5; A4; B1; B2; B3; B4; C1; C2; D2; D3; D4.
- 1 jeringa de 4 g de EPIC Esmalte en los colores: EA1; EA2; EA3; EA3.5; EA4; EB1; EB2; EB3; EB4; EC1; EC2; EC3; EC4; ED2; ED3; ED4.
- 1 jeringa de 4 g de EPIC Dentina en los colores: DA1; DA2; DA3; DA3.5; DA4; DB0.5; DB1; DB2; DB3; DB4; DC1; DC2; DC3; DC4; DD2; DD3; DD4.
- 1 jeringa de 4 g de EPIC Dientes Blanqueados: WB; WO; WE; BL1 y BL2.
- 1 jeringa de 4 g de EPIC Opaquer en el color A1.
- 1 jeringa de 4 g de EPIC Effects: Translúcido e Incisal.
- 1 jeringa de 2 g de EPIC Pigments: Amarillo; Azul; Blanco Marrón; Ocre and Violeta.
- 1 jeringa de 2 g de EPIC Flow en los colores Opalescente y Opaquer B1.

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